

# IP6K1 Antibody



PACO09979

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

IP6K1

**Uniprot**

Q92551

**Synonyms:**

inositol hexakisphosphate kinase 1; IP6K1; IHPK1; MGC9925; PiUS ;

**Immunogen:**

Human IP6K1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A