

PACO09984

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

IPPK

Uniprot

Q9H8X2

Synonyms:

inositol 1,3,4,5,6-pentakisphosphate 2-kinase; IPPK; C9orf12; FLJ13163; INSP5K2; IP5K; IPK1; KIAA0699; bA476B13.1 ;

Immunogen:

Human IPPK.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A