

IRF2BP1 Antibody



PACO09987

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

IRF2BP1

Uniprot

Q8IU81

Synonyms:

interferon regulatory factor 2 binding protein 1; IRF2BP1; DKFZp434M154 ;

Immunogen:

Human IRF2BP1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A