

IRF3 Antibody



PACO09989

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

IRF3

Uniprot

Q14653

Synonyms:

interferon regulatory factor 3;

Immunogen:

Human IRF3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A