

IRGC Antibody



PACO09994

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

IRGC

Uniprot

Q6NXR0

Synonyms:

immunity-related GTPase family, cinema; IRGC; CINEMA; IRGC1; ligp5; R30953_1 ;

Immunogen:

Human IRGC.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A