

ITGB1BP2 Antibody



PACO10008

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

ITGB1BP2

Uniprot

Q9UKP3

Synonyms:

integrin beta 1 binding protein (melusin) 2; ITGB1BP2; CHORDC3; ITGB1BP; MELUSIN; MGC119214; MSTP015 ;

Immunogen:

Human ITGB1BP2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A