

IWS1 Antibody



PACO10023

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

IWS1

Uniprot

Q96ST2

Synonyms:

IWS1 homolog (S. cerevisiae); IWS1; DKFZp761G0123; FLJ10006; FLJ14655; FLJ32319; MGC126375; MGC126376 ;

Immunogen:

Human IWS1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A