

IYD Antibody



PACO10024

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

IYD

Uniprot

Q6PHW0

Synonyms:

iodotyrosine deiodinase; IYD; C6orf71; DEHAL1; dJ422F24.1 ;

Immunogen:

Human IYD.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A