

JAKMIP2 Antibody



PACO10026

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

JAKMIP2

Uniprot

Q96AA8

Synonyms:

janus kinase and microtubule interacting protein 2; JAKMIP2; JAMIP2; KIAA0555; NECC1 ;

Immunogen:

Human JAKMIP2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A