

JAM2 Antibody



PACO10027

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

JAM2

Uniprot

P57087

Synonyms:

junctional adhesion molecule 2; JAM2; C21orf43; CD322; JAM-B; JAMB; PRO245; VE-JAM; VEJAM ;

Immunogen:

Human JAM2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A