

JPH2 Antibody



PACO10029

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

JPH2

Uniprot

Q9BR39

Synonyms:

junctophilin 2; Jph2;1110002E14Rik; JP-2; Jp2 ;

Immunogen:

Human Jph2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A