

# AJUBA Antibody



PACO10031

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

AJUBA

**Uniprot**

Q96IF1

**Synonyms:**

jub, ajuba homolog (Xenopus laevis); JUB; Ajuba; MGC15563 ;

**Immunogen:**

Human JUB.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A