

KALRN Antibody



PACO10034

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

KALRN

Uniprot

O60229

Synonyms:

kalirin, RhoGEF kinase; KALRN; DUET; FLJ16443; HAPIP; TRAD; duo ;

Immunogen:

Human KALRN.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A