

KANK1 Antibody



PACO10035

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

KANK1

Uniprot

Q14678

Synonyms:

KN motif and ankyrin repeat domains 1; KANK1; ANKRD15; DKFZp451G231; KANK; KIAA0172; MGC43128 ;

Immunogen:

Human KANK1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A