

KAT2A Antibody



PACO10038

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

KAT2A

Uniprot

Q92830

Synonyms:

K(lysine) acetyltransferase 2A; KAT2A; GCN5; GCN5L2; MGC102791; PCAF-b; hGCN5 ;

Immunogen:

Human KAT2A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A