

KAT2B Antibody



PACO10039

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

KAT2B

Uniprot

Q92831

Synonyms:

K(lysine) acetyltransferase 2B; KAT2B; CAF; P; P/CAF; PCAF ; P300/CBP associated factor

Immunogen:

Human KAT2B.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A