

KATNB1 Antibody



PACO10042

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

KATNB1

Uniprot

Q9BVA0

Synonyms:

katanin p80 (WD repeat containing) subunit B 1; KATNB1; KAT ;

Immunogen:

Human KATNB1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A