

# KCNAB1 Antibody



PACO10045

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

KCNAB1

**Uniprot**

Q14722

**Synonyms:**

potassium voltage-gated channel, shaker-related subfamily, beta member 1; KCNAB1; AKR6A3; KCNA1B; KV-BETA-1; Kvb1.3; hKvBeta3; hKvb3 ;

**Immunogen:**

Human KCNAB1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A