

KCNF1 Antibody



PACO10051

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

KCNF1

Uniprot

Q9H3M0

Synonyms:

potassium voltage-gated channel, subfamily F, member 1; KCNF1; IK8; KCNF; KV5.1; MGC33316; kH1 ;

Immunogen:

Human KCNF1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A