

KCNIP4 Antibody



PACO10055

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

KCNIP4

Uniprot

Q6PIL6

Synonyms:

Kv channel interacting protein 4; KCNIP4; CALP; KCHIP4; MGC44947 ;

Immunogen:

Human KCNIP4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A