

KCNN1 Antibody



PACO10067

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

KCNN1

Uniprot

Q92952

Synonyms:

potassium intermediate/small conductance calcium-activated channel, subfamily N, member 1; KCNN1; KCa2.1; SK1; SKCA1; hSK1 ;

Immunogen:

Human KCNN1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A