

KCTD14 Antibody



PACO10070

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

KCTD14

Uniprot

Q9BQ13

Synonyms:

potassium channel tetramerisation domain containing 14; KCTD14; MGC2376 ;

Immunogen:

Human KCTD14.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A