

KCTD5 Antibody



PACO10072

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

KCTD5

Uniprot

Q9NXV2

Synonyms:

potassium channel tetramerisation domain containing 5; KCTD5; FLJ20040 ;

Immunogen:

Human KCTD5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A