

PACO10077

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

KDSR

Uniprot

Q06136

Synonyms:

3-ketodihydrosphingosine reductase; KDSR; DHSR; FLJ36555; FLJ92680; FVT1; SDR35C1 ;

Immunogen:

Human KDSR.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A