

KHDRBS1 Antibody



PACO10079

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

KHDRBS1

Uniprot

Q07666

Synonyms:

KH domain containing, RNA binding, signal transduction associated 1; KHDRBS1; FLJ34027; Sam68; p62 ;

Immunogen:

Human KHDRBS1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A