

KIAA0319L Antibody



PACO10082

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

KIAA0319L

Uniprot

Q8IZA0

Synonyms:

KIAA0319-like; KIAA0319L; KIAA1837; PKD1-like; PP791; RP4-765A10.3 ;

Immunogen:

Human KIAA0319L

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A