

KIAA1191 Antibody



PACO10087

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

KIAA1191

Uniprot

Q96A73

Synonyms:

KIAA1191; KIAA1191; FLJ21022; p60MONOX ;

Immunogen:

Human KIAA1191.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A