

KIF19 Antibody



PACO10095

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:**Protein Background:****Gene ID:**

KIF19

Uniprot

Q2TAC6

Synonyms:

kinesin family member 19; KIF19; FLJ37300; KIF19A ;

Immunogen:

Human KIF19.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A