

KIFC1 Antibody



PACO10109

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

KIFC1

Uniprot

Q9BW19

Synonyms:

kinesin family member C1; KIFC1; HSET; KNSL2; MGC1202; MGC149736; MGC149737 ;

Immunogen:

Human KIFC1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A