

KLF13 Antibody



PACO10122

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

KLF13

Uniprot

Q9Y2Y9

Synonyms:

Kruppel-like factor 13; KLF13; BTEB3; FKLF2; NSLP1; RFLAT-1; RFLAT1 ;

Immunogen:

Human KLF13.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A