

KLHDC4 Antibody



PACO10128

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

KLHDC4

Uniprot

Q8TBB5

Synonyms:

kelch domain containing 4; KLHDC4; DKFZp434G0522; FLJ00104 ;

Immunogen:

Human KLHDC4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A