

KLHL12 Antibody



PACO10129

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

KLHL12

Uniprot

Q53G59

Synonyms:

kelch-like 12 (Drosophila); KLHL12; C3IP1; DKIR; FLJ27152 ;

Immunogen:

Human KLHL12.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A