

KLHL14 Antibody



PACO10131

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

KLHL14

Uniprot

Q9P2G3

Synonyms:

kelch-like 14 (Drosophila);

Immunogen:

Human KLHL14.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A