

KLHL18 Antibody



PACO10132

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

KLHL18

Uniprot

O94889

Synonyms:

kelch-like 18 (Drosophila); KLHL18; FLJ13703; FLJ61265; KIAA0795 ;

Immunogen:

Human KLHL18.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A