

KLHL2 Antibody



PACO10133

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

KLHL2

Uniprot

O95198

Synonyms:

kelch-like 2, Mayven (Drosophila); KLHL2; ABP-KELCH; MAV; MAYVEN ;

Immunogen:

Human KLHL2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A