

KLHL20 Antibody



PACO10134

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

KLHL20

Uniprot

Q9Y2M5

Synonyms:

kelch-like 20 (Drosophila); KLHL20; KHLHX; KLEIP; KLHLX; RP3-383J4.3 ;

Immunogen:

Human KLHL20.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A