

KLHL22 Antibody



PACO10136

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

KLHL22

Uniprot

Q53GT1

Synonyms:

kelch-like 22 (Drosophila); KLHL22; KELCHL ;

Immunogen:

Human KLHL22.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A