

KLHL35 Antibody



PACO10140

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

KLHL35

Uniprot

Q6PF15

Synonyms:

kelch-like 35 (Drosophila); KLHL35;26597; FLJ33790 ;

Immunogen:

Human KLHL35.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A