

KLHL4 Antibody



PACO10141

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

KLHL4

Uniprot

Q9C0H6

Synonyms:

kelch-like 4 (Drosophila); KLHL4; DKELCHL; KHL4; KIAA1687 ;

Immunogen:

Human KLHL4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A