

KLRG1 Antibody



PACO10148

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

KLRG1

Uniprot

Q96E93

Synonyms:

killer cell lectin-like receptor subfamily G, member 1; KLRG1;2F1; CLEC15A; MAFA; MAFA-2F1; MAFA-L; MAFA-LIKE; MGC13600 ;

Immunogen:

Human KLRG1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A