

KRT33B Antibody



PACO10168

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

KRT33B

Uniprot

Q14525

Synonyms:

keratin 33B; KRT33B; HA3II; Ha-3II; KRTHA3A; KRTHA3B; hHa3-II ;

Immunogen:

Human KRT33B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A