

# KRT34 Antibody



PACO10169

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

KRT34

**Uniprot**

O76011

**Synonyms:**

436aa,49kd;

**Immunogen:**

Human KRT34.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20oC, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A