

LACTB Antibody



PACO10183

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

LACTB

Uniprot

P83111

Synonyms:

lactamase, beta; LACTB; FLJ14902; G24; MRPL56 ;

Immunogen:

Human LACTB.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A