

LAMC2 Antibody



PACO10188

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

LAMC2

Uniprot

Q13753

Synonyms:

laminin, gamma 2; LAMC2; B2T; BM600; CSF; EBR2; EBR2A; LAMB2T; LAMNB2; MGC138491; MGC141938 ;

Immunogen:

Human LAMC2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A