

LARP7 Antibody



PACO10196

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

LARP7

Uniprot

Q4G0J3

Synonyms:

La ribonucleoprotein domain family, member 7; LARP7; DKFZp564K112; HDCMA18P; MGC104360; PIP7S ;

Immunogen:

Human LARP7.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A