

LARS2 Antibody



PACO10198

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

LARS2

Uniprot

Q15031

Synonyms:

leucyl-tRNA synthetase 2, mitochondrial; LARS2; KIAA0028; LEURS; MGC26121 ;

Immunogen:

Human LARS2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A