

LAT2 Antibody



PACO10202

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

LAT2

Uniprot

Q9GZY6

Synonyms:

linker for activation of T cells family, member 2; LAT2; HSPC046; LAB; NTAL; WBSCR15; WBSCR5; WSCR5 ;

Immunogen:

Human LAT2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A