

LCE1B Antibody



PACO10209

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

LCE1B

Uniprot

Q5T7P3

Synonyms:

late cornified envelope 1B; LCE1B; LEP2; SPRL2A ;

Immunogen:

Human LCE1B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A