

PACO10211

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

LCOR

Uniprot

Q96JN0

Synonyms:

ligand dependent nuclear receptor corepressor; LCOR; FLJ38026; KIAA1795; MLR2; RP11-175O19.1 ;

Immunogen:

Human LCOR.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A