

LDLRAD1 Antibody



PACO10223

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

LDLRAD1

Uniprot

Q5T700

Synonyms:

low density lipoprotein receptor class A domain containing 1;

Immunogen:

Human LDLRAD1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A